

Laser Desorption/Ionization Time-of-Flight Mass Spectrometry for Future *In Situ* Planetary Missions



Stephanie A. Getty, William B. Brinckerhoff, Xiang Li, Melissa A. M. Floyd,
Ricardo D. Arevalo, Jr., Jamie Elsila, and Michael Callahan

NASA Goddard Space Flight Center

Timothy J. Cornish and Scott A. Ecelberger

C&E Research, Inc.

Nancy Chanover, David Voelz, Xifeng Xiao, Rula Tawalbeh, Kyle Uckert, and David Glenar


New Mexico State University

IN SITU INSTRUMENT DEVELOPMENT: LD-TOF-MS IN A COMPACT INSTRUMENT



Commercial LD-TOF-MS is a gold-standard technique for the analysis of biomolecules

We have built a miniature LD-TOF-MS smaller by an order of magnitude




IN SITU INSTRUMENT DEVELOPMENT: LD-TOF-MS FOR IN SITU COMPOSITION


LD-TOF is more generally capable of broadband detection of inorganic and organic composition directly from a solid sample

- Mineral composition
- Small organics: amino acids, carboxylic acids, polycyclic aromatics, etc.
- Intermediate organics: molecular fossil precursors, conjugated polymers, etc.
- Large organics: peptides, biopolymers, etc.
- Can resolve elemental isotopes and ^{12}C , ^{13}C isotope patterns in organics

We have demonstrated this broad detection capability using a diverse suite of standards, analogs, and simulants



INSTRUMENT DESIGN: 5 KG-CLASS IN SITU ANALYZER



LD-TOF-MS Mass Estimate	
Subsystem	Mass/g
1. TOF-MS	700
1.1 Mass Analyzer	330
1.2 Housing	370
2. Laser	1693
3. Optical	380
4. Electronics	1826
4.1 Comm/Data	291
4.2 Power Supply	585
4.3 Pulsed HV	440
4.4 Detector	260
4.5 Harness	250
5. Turbo Pump	550
5.1 Pump, 200 krpm	200
5.2 Controller	350
TOTAL	5149
TOTAL + 30% Margin	6848

COMPARISON: ADVANTAGES OF LD-TOF-MS

	GC-QMS (SAM)	LDI-ITMS (MOMA-MS)	LD-TOF-MS
Sample: Volatile vs Non-volatile	Volatile	Non-volatile (Volatile with GC- EI)	Non-volatile
Power Supply	RF	RF	Pulsed DC
Mass Range	2-535 Da	1-2,000 Da	1-150,000 Da
Ion Polarity	Positive ions	Positive ions	Positive and Negative
Instrument Mass [Estimates]	~15 kg	~7 kg	~5 kg

EXAMPLE SPECTRA:

STANDARDS AND ANALOG SAMPLES

Positive Mode

- Inorganic (cation) Minerals •
 - Polycyclic Aromatics •
- High Molecular Weight Polymers and Biopolymers •

Negative Mode

- Inorganic (anion) minerals •
- Amino Acids • Carboxylic Acids •
 - Analog Mixtures •

Integrated *In Situ* Instrumentation

- Correlated IR and Mass Spectra •
- Advantages in Surface Operations, Confidence in Analyses •

Advanced Mass Spectrometry

- Molecular Fingerprinting through Fragment Analysis •
- Specificity through Two-Step Laser Desorption/Ionization (L2MS) •

EXAMPLE SPECTRA:

STANDARDS AND ANALOG SAMPLES

Positive Mode

- Inorganic (cation) Minerals •
 - Polycyclic Aromatics •
- High Molecular Weight Polymers and Biopolymers •

Negative Mode

- Inorganic (anion) minerals •
- Amino Acids • Carboxylic Acids •
- Analog Mixtures •

Integrated *In Situ* Instrumentation

- Correlated IR and Mass Spectra •
- Advantages in Surface Operations, Confidence in Analyses •

Advanced Mass Spectrometry


- Molecular Fingerprinting through Fragment Analysis •
- Specificity through Two-Step Laser Desorption/Ionization (L2MS) •

POSITIVE ION MODE:

INORGANIC (CATION) COMPOSITION

Sedimentary and Aqueously Altered Minerals,
e.g., Sulfates, Carbonates, Iron Oxides, etc.

MCP Signal (V)




COMPLEMENTARY POSITIVE AND NEGATIVE ION DETECTION: E.G., PERCHLORATES AND SULFIDES

Sodium Perchlorate (NaClO_4)


Positive Ion Mode

(+) Positive Ion



Negative Ion Mode


(-) Negative Ion




Pyrite (FeS_2)

Positive Ion Mode


Pyrite: FeS_2



Negative Ion Mode




COMPLEMENTARY POSITIVE AND NEGATIVE ION DETECTION: DETECTING ORGANICS ACROSS CLASSES




ORGANIC COMPOSITION

Beta-Carotene: relevant compound to molecular fossil precursors,
polyacetylene synthesis in interstellar clouds



EXAMPLE MEASUREMENTS:

BIOMOLECULES NEXT GENERATION *IN SITU* ASTROBIOLOGY



EXAMPLE SPECTRA:

STANDARDS AND ANALOG SAMPLES

Positive Mode

- Inorganic (cation) Minerals •
 - Polycyclic Aromatics •
- High Molecular Weight Polymers and Biopolymers •

Negative Mode

- Inorganic (anion) minerals •
- Amino Acids • Carboxylic Acids •
- Analog Mixtures •

Integrated In Situ Instrumentation







- Correlated IR and Mass Spectra •
- Advantages in Surface Operations, Confidence in Analyses •

Advanced Mass Spectrometry

- Molecular Fingerprinting through Fragment Analysis •
- Specificity through Two-Step Laser Desorption/Ionization (L2MS) •

COMPLEMENTARY IN SITU ANALYSES: COINCIDENT IR AND MASS SPECTRA

NMSU IR Spectrometer



EXAMPLE SPECTRA:

STANDARDS AND ANALOG SAMPLES

Positive Mode

- Inorganic (cation) Minerals •
 - Polycyclic Aromatics •
- High Molecular Weight Polymers and Biopolymers •

Negative Mode

- Inorganic (anion) minerals •
- Amino Acids • Carboxylic Acids •
- Analog Mixtures •

Integrated *In Situ* Instrumentation


- Correlated IR and Mass Spectra •
- Advantages in Surface Operations, Confidence in Analyses •

Advanced Mass Spectrometry

- Molecular Fingerprinting through Fragment Analysis •
- Specificity through Two-Step Laser Desorption/Ionization (L2MS) •


PRECISION ION GATING AND TANDEM MS

Structural determination using MS/MS techniques



ADVANCED MASS SPECTROMETRY: TWO-STEP LASER IONIZATION FOR SPECIES SELECTIVITY IN A COMPLEX MIXTURE

Isolate the non-polar, aromatic fraction of complex sample
without front-end chromatography




LD-TOF-MS offers a compelling *in situ* capability for exploration of high priority planetary surface materials: Sediments – Regolith – Ices
Contextual Mineralogy • Organic Building Blocks • Biological Polymers

Planetary Evolution:

Evidence for extinct/extant life
Preserved organics
Sample selection
Mineralogical context

**Sediments
and Salts**



Primitive Materials:

Prebiotic chemistry
Organics available to early Earth
Sample selection
Mineralogical context

